



Internal Watercooler Sanitiser

Product codes & Pack size: AQC004 (5litre) AQC003 (1litre)

Where to use: To sanitise the internal water tanks and components on water coolers, product is dosed into machine via inlet water piper or poured directly into the reservoir tank

Key Benefits v Competition

- Unique colour indicator tells customer when product has flowed through the system and visually confirms all parts are being sanitised
- Colour indicator also confirms when product has been flushed out of the watercooler, make it safe and easy to use. No need for test strips.
- 3 stage all in one mild descaling, cleans and sanitising formula.
- Effective against bio-film as well as bacteria.
- Cleaning card available to recommend cleaning regimes.

Usage/Dilution rates:

Should be used as part of regular sanitising and cleaning programme every 3-4 months. 25ml of solution for every 2 litres, typical tanks contain 2 litres of water but need to check watercooler tank spec in case larger tank used.



Sales or technical questions

Contact: Nigel Barron internally or Stuart Hay at Bioguard
stuarthay@bioguardhygiene.co.uk 07887 850146